## UNIT-DOSE D-A - MEDICATION TROLLEY COD. 6053/DA

BODY Structure made of powder-coated steel RAL 7035, thickness 10/10, finished with acrylic polyurethane enamels resistant to scratches, stains and corrosion. The structure can be washed and disinfected.

## WORKTOP

Worktop made of self-extinguishing ABS with raised edges on three sides. The worktop is resistant to scratches and stains and can be cleaned, washed and disinfected.

## DRAWER UNIT

Upper part: 8 drawers designed for the positioning of medication bins for unit-dose therapy. Drawers are made of moulded technopolymer with integrated handle and run on soft-closing telescopic slides allowing their complete extraction.

Standard configuration: 32 medication bins (dim. 94x400x82 h mm). Each medication bin is provided with name label holder and adjustable dividers.

Lower part: 1 drawer made of powder-coated steel provided with ergonomic and anti-trauma aluminium handle and stoppers. Standard configuration:1x 200 h mm. Centralized key lock.

#### LEFT SIDE

Ergonomic and anti-trauma aluminium push handle with plastic end caps.

## BASE

Double shock absorbent bumper to prevent the trolley from turning over. 4 thermoformed swivel castors, Ø 125 mm, mounted on bearings; 2 are lockable.

\*Accessories for this trolley are available upon request.\*

Dim.610x600x1300 h mm

Position 1



green blue yellow red white

Rev. 01 del 28/12/2022

# TECHNICAL SHEETS - NR. / -